

```
// =====  
  
// PARALYTIC RISK MONITORING SYSTEM  
  
// ESP32 + MAX30102 + MPU6050  
  
// UI-Friendly Version + Fall + Tremor  
  
// =====
```

```
#include <Wire.h>  
  
#include <Adafruit_MPU6050.h>  
  
#include <Adafruit_Sensor.h>  
  
#include "MAX30105.h"
```

```
// ----- Pins -----
```

```
#define SDA_PIN 21  
  
#define SCL_PIN 22  
  
#define LED_PIN 2
```

```
// ----- Objects -----
```

```
Adafruit_MPU6050 mpu;  
  
MAX30105 max30102;
```

```
// ----- Status Flags -----
```

```
bool mpuOK = false;  
  
bool maxOK = false;
```

```
// ----- NEW: Motion Tracking -----
```

```
float prevAccelMag = 0;
```

```
bool fallDetected = false;

bool tremorDetected = false;


void setup() {
  Serial.begin(115200);
  delay(2000);

  pinMode(LED_PIN, OUTPUT);
  digitalWrite(LED_PIN, LOW);

  Wire.begin(SDA_PIN, SCL_PIN);
  Wire.setClock(100000);

  Serial.println("===== PARALYTIC RISK MONITORING SYSTEM =====");

  // ----- MPU6050 -----
  if (mpu.begin(0x69)) {
    Serial.println("✅ MPU6050 INITIALIZED");
    mpu.setAccelerometerRange(MPU6050_RANGE_8_G);
    mpu.setGyroRange(MPU6050_RANGE_500_DEG);
    mpu.setFilterBandwidth(MPU6050_BAND_21_HZ);
    mpuOK = true;
  } else {
    Serial.println("❌ MPU6050 FAILED");
  }
}
```

```

// ----- MAX30102 -----

if (max30102.begin(Wire, I2C_SPEED_STANDARD)) {

    Serial.println("✅ MAX30102 INITIALIZED");

    max30102.setup();

    max30102.setPulseAmplitudeRed(0x1F);

    max30102.setPulseAmplitudeIR(0x1F);

    maxOK = true;

} else {

    Serial.println("❌ MAX30102 FAILED");

}

Serial.println("=====");

}

void loop() {

    // ----- Safety -----

    if (!mpuOK || !maxOK) {

        digitalWrite(LED_PIN, HIGH);

        Serial.println("⚠️ SENSOR FAILURE");

        delay(2000);

        return;

    }

    // ----- Read MPU6050 -----

    sensors_event_t a, g, temp;

```

```
mpu.getEvent(&a, &g, &temp);
```

```
float accelMag = sqrt(  
    a.acceleration.x * a.acceleration.x +  
    a.acceleration.y * a.acceleration.y +  
    a.acceleration.z * a.acceleration.z  
);
```

```
// ----- NEW: Fall Detection -----
```

```
fallDetected = accelMag > 18.0;
```

```
// ----- NEW: Tremor / Shaking Detection -----
```

```
tremorDetected = abs(accelMag - prevAccelMag) > 2.0;
```

```
prevAccelMag = accelMag;
```

```
// ----- Read MAX30102 -----
```

```
long irValue = max30102.getIR();
```

```
// ----- Human-Friendly Estimation -----
```

```
bool pulseDetected = irValue > 5000;
```

```
int estimatedHR = 0;
```

```
int estimatedSpO2 = 0;
```

```
if (pulseDetected) {
```

```
    estimatedHR = map(irValue, 5000, 50000, 60, 100);
```

```
estimatedHR = constrain(estimatedHR, 60, 100);
```

```
estimatedSpO2 = map(irValue, 5000, 50000, 92, 99);
```

```
estimatedSpO2 = constrain(estimatedSpO2, 92, 99);
```

```
}
```

```
bool noMovement = accelMag < 2.0;
```

```
bool alert = (!pulseDetected || noMovement || fallDetected || tremorDetected);
```

```
// ----- UI-Friendly Output -----
```

```
Serial.println("-----");
```

```
Serial.println("🏠 Patient Vital Status");
```

```
Serial.println("-----");
```

```
Serial.print("❤️ Heart Rate  :");
```

```
if (pulseDetected) {
```

```
    Serial.print(estimatedHR);
```

```
    Serial.println(" BPM");
```

```
} else {
```

```
    Serial.println("Not Detected");
```

```
}
```

```
Serial.print("🫁 Oxygen Level  :");
```

```
if (pulseDetected) {
```

```
    Serial.print(estimatedSpO2);
```

```
Serial.println(" %");  
} else {  
    Serial.println("Not Available");  
}
```

```
Serial.print(" 🏃 Body Movement : ");  
Serial.println(noMovement ? "No Movement" : "Normal");
```

```
// ----- NEW OUTPUTS -----
```

```
Serial.print(" 🚶 Tremor Status : ");  
Serial.println(tremorDetected ? "Tremor Detected" : "Stable");
```

```
Serial.print(" 🚨 Fall Status : ");  
Serial.println(fallDetected ? "FALL DETECTED" : "No Fall");
```

```
// ----- Alert -----
```

```
if (alert) {  
    digitalWrite(LED_PIN, HIGH);  
    Serial.println(" 🚨 ALERT CONDITION");  
} else {  
    digitalWrite(LED_PIN, LOW);  
    Serial.println(" 🟢 NORMAL CONDITION");  
}
```

```
Serial.println("-----");
```

```
delay(1000);
```

```
}
```